






Arthrosmart Anchor[®]





Features

- Titanium self-tapping screw (6Al4V alloy), providing increased pullout resistance;
- Disposable applicator, facilitating the insertion technique;
- Mounted with two non-absorbable suture strands (Fiber®), with high resistance rupture;
- Screw-type anchors available in diameters of 2.0 mm, 2.7 mm, 3.5 mm, 4.0 mm and 5.0 mm for indication of use in ATM, hand, wrist, elbow, shoulder, foot, ankle, knee and hip.

Indications of Use

	ATM	Hand / Elbow	Hand	Shoulder	Foot / Ankle	Knee	Hip
Ø 2,0	X	X	-	-	-	-	-
Ø 2,7	X	X	X	X	X	X	X
Ø 3,5	-	X	X	X	X	X	X
Ø 4,0	-	-	-	X	X	X	X
Ø 5,0	-	-	-	X	X	X	X

Information extracted from the Product Use Instruction

Anhcor with Disposable Applicator



TiAlV

Code

39.44.04.00020
39.44.04.00027
39.44.04.00035
39.44.04.00040
39.44.04.00050

Ø Anchor

Ø 2,0 mm
Ø 2,7 mm
Ø 3,5 mm
Ø 4,0 mm
Ø 5,0 mm

Applicator

Standard
Standard
Standard
Standard
Standard



with Needle

39.44.03.10020
39.44.03.10027
39.44.03.00035
39.44.03.00040
39.44.03.00050

Ø 2,0 mm
Ø 2,7 mm
Ø 3,5 mm
Ø 4,0 mm
Ø 5,0 mm

Short
Short
Standard
Standard
Standard